

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS

ASMS 2018 WP-249

Sven Vedder¹, Anja Grüning¹, Julia Sander¹
1 Shimadzu Europa, Albert-Hahn-Str. 6-10,
47269 Duisburg, Germany

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS

Introduction

Pyrrolizidine alkaloids (PAs) are secondary plant metabolites that are supposed to be carcinogenic and genotoxic. They occur mainly in plants of the Boraginaceae, Asteraceae and Fabaceae families. They contain a pyrrolizidine core and make up a large group of heterocyclic alkaloids mainly derived from the 4 Necin bases platynecine, retronecine, heliotridin and ontonecin. PAs are hepatotoxic if they carry a 1,2-double bond as well as an esterified side chain which is a structural prerequisite for their hepatic activation. Plant food and beverage, phytopharmaceuticals or even

animal feed can easily be contaminated with PAs and enter the food chain. Currently there are discussions on possible regulatory measures caused by the presence of PAs in honey, tea, herbal infusions and food supplements. Existing methods include laborious sample preparation, e.g. solid-liquid extraction followed by solid phase extraction for clean-up. Here we report an on-line SPE UHPLC-MS/MS method, which overcomes the difficulties of combining low pressure online SPE with high pressure analytical UHPLC.




Figure 1 General structure of pyrrolizidine alkaloids

Methods and Materials

Sample Preparation

Tea samples were extracted twice with 0.05M sulfuric acid by sonication. Before centrifugation the pH of the combined extracts was adjusted with ammonium hydroxide.

UHPLC method

| | |
|--------------------|--|
| Instrument | : Nexera UHPLC, Shimadzu |
| Column | : Shim-pack XR-ODS III, 150 mm x 2.0 mm, 2.2 μm, Shimadzu |
| Mobile phase A | : 5 mM ammonium formate + 0.1% formic acid |
| B | : methanol + 5 mM ammonium formate + 0.1% formic acid |
| Flow rate | : 0.4 mL/min |
| Time program | : B conc. 1% (0-1.6 min) -50% (14.6 min) – 71.5% (18.1 min) – 95% (18.2 min – 20.2 min) -1% (20.3 min – 25 min) |
| Column temperature | : 30 °C |

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS

Online SPE method

| | |
|--------------------|--|
| Column | : EVOLUTE® EXPRESS ABN, 30 x 2.1 mm, Biotage |
| Mobile phase | : 5 mM ammonium formate + 0.1% formic acid for sample loading methanol / H ₂ O + 5 mM ammonium formate + 0.1% formic acid methanol, isopropanol for washing of SPE column |
| Flow rate | : 0.2 / 2 mL/min |
| Injection vol. | : 50 µL |
| Column temperature | : RT |

MS conditions

| | |
|------------------------|-----------------------|
| Instrument | : LCMS-8060, Shimadzu |
| Ionization | : pos ESI |
| Nebulizing gas | : 3 L/min |
| Heating gas | : 15 L/min |
| Drying gas | : 5 L/min |
| Interface temperature | : 400 °C |
| DL temperature | : 300 °C |
| Heat block temperature | : 400 °C |
| CID gas | : 270 kPa |
| Interface voltage | : 1 kV |

Result

Method development of the online SPE

The neutralized and centrifuged tea extract samples were put into the autosampler and transferred to the on-line SPE column using an aqueous solution. After washing the sample was eluted with only 10 µL solvent and trapped into a loop. By switching the loop the eluted sample was transferred to the analytical column. A binary gradient

separated the PAs for quantification. Due to this hardware set-up UHPLC with high backpressure and on-line SPE which is pressure limited were successfully combined. By careful fine-tuning of the SPE elution and the chromatographic conditions the separation of critical peak pairs could be maintained.




Figure 2 Typical chromatogram of pyrrolizidine alkaloids in tea matrix including the separated pressure curves of the analytical column (Pump A and B pressure) and the online SPE column (Pump C pressure)

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS




Figure 3 Setup of the on-line SPE analytical system

Quantitative Analysis of tea samples

By using the reported instrument set-up, analysis and thus the quantification of 16 PAs and 14 of their related N-Oxides could be performed. Calibration curves in different tea matrices (black tea, green tea and herbal tea) determined in duplicate showed good precision and accuracy and even in a complex matrix like tea we were able to easily quantify the PAs in at least the range of 10 to 400 µg/kg. This is comparable to the established methods

using manual sample preparation. For all analytes, weighted regression resulting in r^2 0.99 could be achieved, with $S/N > 10$ for LLOQ levels.

Exemplary calibration curves obtained for the 30 compounds are shown in Figure 4, Chromatograms of exemplary LLOQs are shown in Figure 5, the LLOQs which could be achieved in the different tea matrices are shown in Table 1.

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS




Figure 4 Exemplary calibration curves in different tea matrices




Figure 5 Exemplary chromatograms of LLOQs different tea matrices

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS

Table 1 LLOQs of the pyrrolizidine alkaloids in different tea matrices

| | Black tea LLOQ | | Green tea LLOQ | | Herbal tea LLOQ | |
|--|----------------|-------|----------------|-------|-----------------|-------|
| | ng/mL | µg/kg | ng/mL | µg/kg | ng/mL | µg/kg |
| Echimidine | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Echimidine-N-oxide | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Erucifoline | 0.20 | 4.0 | 0.20 | 4.0 | 0.50 | 10.0 |
| Erucifoline-N-Oxide | 0.10 | 2.0 | 0.20 | 4.0 | 0.10 | 2.0 |
| Europine | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Europine -N-Oxide | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Heliotrine | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Heliotrine N-oxide | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Intermedine | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Intermedine N-oxide / Indicine N-oxide | N.A | N.A | 0.05 | 1.0 | 0.05 | 1.0 |
| Jacobine | 0.20 | 4.0 | 0.50 | 10.0 | 0.50 | 10.0 |
| Jacobine N-oxide | 0.05 | 1.0 | 0.10 | 2.0 | 0.10 | 2.0 |
| Lasiocarpine | 0.05 | 1.0 | 0.10 | 2.0 | 0.05 | 1.0 |
| Lasiocarpine N-oxide | 0.10 | 2.0 | 0.10 | 2.0 | 0.05 | 1.0 |
| Lycopsamine / Indicine | N.A | N.A | 0.05 | 1.0 | 0.05 | 1.0 |
| Lycopsamine N-oxid | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Monocrotaline | 0.50 | 10.0 | 0.20 | 4.0 | 0.50 | 10.0 |
| Monocrotaline-N-oxide | 0.10 | 2.0 | 0.10 | 2.0 | 0.10 | 2.0 |
| Retrosine | 0.10 | 2.0 | 0.05 | 1.0 | 0.10 | 2.0 |
| Retrosine N-oxide | 0.20 | 4.0 | 0.50 | 10.0 | 0.50 | 10.0 |
| Senecionine | 0.05 | 1.0 | 0.50 | 10.0 | 0.50 | 10.0 |
| Senecionine N-oxide | 0.05 | 1.0 | 0.20 | 4.0 | 0.50 | 10.0 |
| Senecyphilline | 0.10 | 2.0 | 0.20 | 4.0 | 0.50 | 10.0 |
| Senecyphilline N-oxide | 0.05 | 1.0 | 0.20 | 4.0 | 0.10 | 2.0 |
| Senecivernine | 0.20 | 4.0 | 0.50 | 10.0 | 0.20 | 4.0 |
| Senecivernine N-oxide | 0.05 | 1.0 | 0.20 | 4.0 | 0.10 | 2.0 |
| Senkirkine | 0.05 | 1.0 | 0.05 | 1.0 | 0.05 | 1.0 |
| Trichodesmine | 0.05 | 1.0 | 0.10 | 2.0 | 0.20 | 4.0 |

A total of 29 commercially available tea samples were analyzed. Among these samples there were 6 samples of green tea, 10 samples of black tea and 13 samples of herbal tea. In 59% of all analyzed tea samples one or more

of the pyrrolizidine alkaloids could be detected above their LLOQ. 3 out of 6 green tea samples, 5 out of 10 black tea samples and 9 out of 13 herbal tea samples were contaminated with pyrrolizidine alkaloids.

Determination of pyrrolizidine alkaloids in plant material using on-line SPE coupled to UHPLC-MS/MS

Conclusions

An on-line SPE method for high-sensitivity analysis was successfully developed for PA analysis in plant material. The manual sample preparation could be reduced to a minimum as the set up of on-line SPE followed by UHPLC-MS/MS saves additional clean-up steps without compromising the performance of the assay.

The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures

First Edition: June, 2018

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.